


Characterization of microporous carbon anodes using N₂ and CO₂ adsorption isotherms on TriStar II Plus 3030



TRISTAR

The rapid evolution of silicon-carbon (Si-C) composite anodes and sodium-ion batteries has placed a spotlight on the importance of microporous carbons. In Si-C anodes, high microporosity provides a robust scaffold for silicon deposition, maximizing specific surface area and capacity. For sodium-ion batteries, hard carbon remains the leading candidate for efficient ion intercalation.

A wide range of carbon precursors can be used to produce porous carbon-based anodes, including coal-derived carbon, pitch-derived carbon, polysaccharide-derived carbon, and biomass-derived carbon. Among these, utilizing biomass waste, such as fruit shells, peels, straw, and rice husk, through high-temperature carbonization and structural restructuring to prepare high-quality microporous anode materials is an economical and environmentally friendly approach that aligns with global resource recycling and utilization policy.

Advanced Micropore Characterization via Dual-Gas Adsorption

The use of microporous carbon anode materials necessitates reliable characterization of micropores using advanced gas adsorption techniques. Standard N₂ gas adsorption instruments often fall short when characterizing the finest micropores due to the absence of low-pressure transducers. To address this, **Micromeritics** utilizes a dual-probe approach on its industry-leading **TriStar II Plus 3030 instruments**:

1. **Nitrogen (N₂) Adsorption:** Provides a broad view of mesopores and larger micropores.
2. **Carbon Dioxide (CO₂) Adsorption** (at 273 K): CO₂ is a more linear and smaller molecule (3.3 Å) compared to N₂ (3.64 Å), enabling its diffusion into smaller pores. At 273K, due to the faster kinetics and the high saturation vapor pressure, CO₂ diffuses more rapidly into narrow carbon pores, probing smaller pore sizes more efficiently.

The use of CO₂ as a probe molecule, together with the advanced NLDFT (Non Local Density Functional Theory) model, enables precise determination of pore size distributions across a broad pore range.

Case Study: Investigation of biomass-derived hard carbon anode materials

N₂ Adsorption isotherm of a biomass- derived hard carbon anode material using a Tristar II Plus 3030 is shown in **Figure 1** below.




Figure 1: N₂ adsorption isotherm of biomass-derived microporous hard carbon anode

This is a typical type Ib adsorption isotherm, indicating the presence of micropores and 2-3 nm mesopores in this material. The overall adsorption capacity is also relatively large, and the material has a high pore surface area. During Bulk testing (**Figure 2**), an even more differentiating result was generated by the Permeability test. Metal Powder C generates a significantly higher Pressure Drop across the powder bed than the other samples, indicating that Metal Powder C is considerably less permeable than Powders A and B.




Figure 2: pore size distribution calculated from N₂ isotherm in biomass-derived microporous hard carbon anode

The pore size distribution calculated using the NDLFT analysis of the N₂ isotherm is shown in **Figure 2**. It reveals pores ranging from 9 Å to 28 Å. Pores below 9 Å are not revealed from this N₂ isotherm. On a TriStar instrument, microporosity can be investigated by collecting an isotherm with CO₂ in addition to the N₂ isotherm.

Tristar allows easy switching up to 3 gases, in addition to a dedicated port for Helium. The CO₂ adsorption isotherm measured at 273 K is shown in **Figure 3**.




Figure 3: CO₂ adsorption isotherm of biomass -derived microporous hard carbon anode

It is well reported in literature that a dual gas NLDFT model can provide a comprehensive analysis of pores in microporous materials¹⁻³. By applying the dual gas NLDFT Advanced PSD method in the MicroActive software – HS-2D-NLDFT for N₂ and NLDFT for CO₂ - we can obtain the complete pore-size distribution covering the micro and mesopore ranges, as shown in **Figure 4**.




Figure 4: Combined pore size distribution of N₂ and CO₂ in biomass-derived microporous hard carbon anode

In this figure, the green curve represents the pore size distribution, with the size corresponding to the peak being called the modal pore size. The blue curve represents the cumulative pore volume. This analysis shows that this material is mainly composed of micropores smaller than 20 Å, and some small mesopores in the 20-30 Å range. This analysis also reveals that the cumulative specific surface area in the pores of this material reaches 1525 m²/g. In microporous materials, the surface area calculated with NLDFT is known to be more accurate than a simple BET estimation.

Porosity Distribution by 2D-NLDFT
 Model: HS-2D-NLDFT, Carbon, N2, 77
 Dual sample: YF39-AC-CVD3-CO2
 Dual model: CO2@273-Carbon, NLDFT
 Isotherm branch: Adsorption
 Method: Non-negative Regularization: 0.03160
 Standard Deviation of Fit: 0.085174 mmol/g

Total Volume in Pores	<=	505.97 Å	:	0.77457 cm ³ /g
Total Area in Pores	>=	3.23 Å	:	1,524.848 m ² /g

Total pore volume and pore surface area of biomass-derived microporous hard carbon anode

Anode material with such a large surface area is a good template to prepare high-silicon-content silicon-carbon anodes, thereby greatly improving the anode's specific capacity.

Model Accuracy and Agreement

As shown in Figure 5, the NLDFT Advanced PSD model shows excellent agreement between the experimental data points and the theoretical adsorption isotherms (solid lines). This validation confirms the reliability of the derived pore volumes and surface area measurements.




Figure 5: Degree of agreement between experimental isotherms and theoretical adsorption isotherms

Conclusion

The biomass-derived hard carbon used in this case study demonstrates exceptionally high surface area, making it an ideal precursor for high-silicon-content anodes. By converting waste into high-value energy storage materials, manufacturers can significantly reduce costs while enhancing specific capacity.

Micromeritics continues to lead in battery material characterization, providing the advanced hardware and analytical models, like the TriStar II Plus 3030 instruments with multiple gas options, and the HS-2D-NLDFT-based analysis model, to reliably analyse the surface area, pore volume, and pore size distribution in these complex microporous materials, unlocking their potential in high-performance energy storage applications.

References

1. <https://pubs.rsc.org/en/content/articlelanding/2022/ma/d2ma00141a>
2. <https://www.sciencedirect.com/science/article/abs/pii/S0008622315004005>
3. <https://www.sciencedirect.com/science/article/abs/pii/S000862231831159X>